ENGLISH

SYMBOLS

MD	Medical Device
MR	Magnetic Resonance (MR) Conditional

INTENDED USE

The device is intended to immobilize the knee. The device must be fitted and adjusted by a healthcare professional.

Indications for use

Situations requiring knee immobilization, for example, traumatic knee injuries, emergency immobilization, and/or pre-/post-op knee immobilization. No known contraindications.

Warnings and Cautions:

- Inspect the skin under the device’s contact area before and after each use.
- The device should be used directly on intact skin and never over an open wound.
- Regular healthcare professional supervision is recommended for patients with peripheral vascular disease, neuropathy, and sensitive skin.
- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.

Non-clinical testing demonstrated that the device is **MR Conditional**. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5 Tesla or 3 Tesla, only.
- Maximum spatial gradient magnetic field of 2,000 Gauss/cm (20 T/m).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
 - If any pain, pain irritation, or unusual reaction occurs with the use of the device.
- The device is for single patient – multiple use.

FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text **(Fig. 1)**.

Device Application

- Unfasten all Straps (A) **(Fig. 2)**.
- Remove the Popliteal Pad (B) and the Kneecap Cover (C).
- If needed, the Posterior Stay (D) may be bent to accommodate the flexion angle of the leg:
 - Place the brace on a table top with one part protruding off the edge. Apply downward pressure to bend the Posterior Stay.
- Verify the correct placement of the Popliteal Pad on the Posterior Stay, and then attach the Popliteal Pad using the double-sided tape on its underside **(Fig. 3)**.
- Carefully lift the patient’s leg and place the device underneath. Align the Posterior Stay along the back of the leg.
- Ensure the kneecap is in line with the notch at the center of the brace. Close the device so that it surrounds the knee, leaving the kneecap exposed **(Fig. 4)**.
- If needed, the Side Panels (E) may be contoured to improve fit and help prevent migration:
 - Place one Side Panel on a table top with one part protruding off the edge. Apply downward pressure to bend the Side Panel to fit the shape of the patient’s leg.
 - Repeat this procedure for the other Side Panel.

- Position and secure both Side Panels to the desired thigh and calf circumference using the hook and loop fastening system **(Fig. 5)**. Once the device has been correctly adjusted, wrap the Side Panels around the leg.
- Secure all six Straps in an alternating above and below fashion until all are secure **(Fig. 6: a, b, c, d, e, f)**.
- Starting below the knee, pull the Strap through the D-ring and secure with the hook fastener (a).
- Re-adjust all Straps for a snug fit, do not overtighten.
- Attach the Kneecap Cover over the knee, if desired **(Fig. 7)**.

Device Removal

- Ask patient to sit down.
- Undo all straps to remove the device from the leg.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and Care

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry, do not expose to direct heat.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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